

Epi-Assist Australia Pty Ltd

ABN 72 075 018 601

SPECIALISING IN EPILEPSY MONITORING & BED EXIT

Reliable Bed Sensor Monitoring Epilepsy Seizures


Reliably monitoring of convulsive seizures is essential for the safety of our loved ones when an epilepsy seizure is taking place.

The Emfit Epileptic Sensor Mat is a state-of-the-art bed monitoring system that will monitor a sleeping adult or child's convulsive seizures from under the mattress.

Epi-Assist Australia Pty Ltd is the sole distributor of the Emfit Safe-Bed, manufactured in Finland. The Emfit Safe-Bed pioneered the Quasi-Piezoelectric sensor, a state-of-the-art bed epileptic and bed exit monitoring system and consists of the L-4060SL bed sensor and the D-1090-2G control unit.

Bed Sensor L-4060SL TGA Certificate: 137139

The L-4060SL patented dynamic thin film ferroelectric sensor mat comes placed between two layers of foam for protection against damage from the bed base and is placed under the mattress. The bed sensor complies with TGA certification standards and has no embedded wires or power running through it. The specially formulated plastic sheet sensor detects a person's seizure in bed by sensing movements and vital body life signs such as heartbeat and breathing.



Bed Sensor L-4060S

HAS NO WEIGHT LIMITATIONS

Control Unit D-1090-2G

The control unit features an ON/OFF/Reset switch. When a seizure is detected an audible alarm with an adjustable volume control can be heard. The alarm tone may be extended to any area of a large home by installing wireless chimes.



Control Unit
D-1090-2G

In the case of adults who wish to retain their independency and remove any burden from relatives, the control unit has a relay output that triggers our medical 'Blue Phone'. Within minutes of a seizure taking place a voice message is relayed to either a landline or mobile telephone.

Emfit offers a manufacturers 2-year limited warranty for both sensors and alarm units.

For further information contact: Epi- Assist Australia Pty Ltd
Phone: Free Call 1800 684 422
Email: sales@modcoms.net
Website: www.emfit.com

All Correspondence to:
P O Box 206
Hunter Regional Mail Centre NSW 2310